

# First Fit card for Nx




# signia

Life sounds brilliant.



The enclosed hearing aids have been First Fit to the Universal Program.

Simply click the icon of the BTE with the green arrow in the tool bar to “Detect connected hearing systems” in Connexx®.



Next, click the radio button “Read out connected instruments” to view the factory-programmed settings in First Fit.

## Next Step



- New Fitting ! All data in the connected instruments, including the audiogram, will be discarded!
- Read out connected instruments
- Mimic fit - transfer session data to connected instruments

**NOTE: Do not click the radio button called “New fitting” or the first fit setting will be cleared. The First Fit Program settings have been optimized to make patient fittings easier and quicker.**

# First Fit information for Nx hearing aid orders\*



## Conduct OVP training with your patient.

Click on the Own Voice Processing tab on the left navigation bar below First Fit.

First Fit



**Own Voice Processing**

Basic Tuning

Prepare Training

Click on "Prepare Training". Connexx provides tips for best OVP results before training starts.

Instruct your patient that when they hear a beep to count continuously aloud in a slightly raised voice for the duration of the training. i.e. 2100, 2200, 2300... This training process usually takes less than 10 seconds.

Start Training

Click on "Start Training".

Have your patient count until the green checkmark appears, indicating that training has been completed.



**Training was successful**

Once training is complete, the "Activate Own Voice Processing" box is automatically checked and the default setting activated in the Settings panel. The default OVP setting is ideal for the vast majority of wearers and should not be adjusted unless it is specifically necessary. Besides the music programs, OVP is active in all other hearing aid programs. If the patient reports that their voice sounds muffled or too soft, consider reducing the setting to min. This may happen for more experienced hearing aid wearers who have become used to hearing their own voice amplified. On the other hand, if a new patient still reports that their voice is too loud, consider adjusting the setting to max.



Activate Own Voice Processing

\*if a *Bluetooth* accessory is ordered the hearing aids have been paired.

The *Bluetooth*\* word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Copyright © 2018 Signia GmbH. All rights reserved. 10240569 5/18 6.0 22B771 11K SI/18525-18