







































Standard BTE Earmolds

												
		Full Shell	Half Shell	Skeleton	Semi-Skeleton ¹	3/4 Skeleton ¹	1/2 Skeleton ¹	Canal	Canal Lock	Open ¹	Open Lock ¹	Skeleton Open ¹
ROES Codes		EM	EM2	EM4	EM5	EM7	EM8	EM9	EM3	EM6	EM10	EM23
Material	Hard	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Soft	✓	✓	✓	-	-	-	✓	✓	-	-	-
Venting	Largest Vent Possible	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Open Air Vent ¹	-	-	-	-	-	-	-	✓	-	-	-
	Select A Vent 4.0mm	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Select A Vent 3.0mm	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Select A Vent 2.0mm	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Semi Iros	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Small Vent	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Power Vent ¹ 1.2mm	✓	✓	✓	✓	✓	✓	-	✓	-	-	-
	Pressure Vent 1.0mm	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	No Vent	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Double Bend Clear 3.0mm OD	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Double Bend Clear 3.3mm OD	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Double Bend Clear 3.5mm OD	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
Tubing	3.1mm and/or Quick-Connect	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	3.3mm and/or Quick-Connect	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Dry 3.3mm and/or Quick-Connect Dry	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	3.6mm and/or Quick-Connect	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Libby Horn	✓	✓	✓	✓	✓	✓	✓	✓	-	-	✓
	Slim Tube Option	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	ThinTube Option	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

	Clear (gloss)		Clear (matte)
	Rose (gloss)		Rose (matte)
	Beige (gloss)		Beige (matte)
	Brown (gloss)		Brown (matte)

	Red Opaque		Blue Opaque		Green Translucent		Red Glitter		Blue Glitter
	Orange Opaque		Lilac Opaque		Pink Translucent		Orange Glitter		Purple Glitter
	Yellow Opaque		Black Opaque		Blue Translucent		Yellow Glitter		Pink Glitter
			White Opaque		Purple Translucent		Green Glitter		Clear Glitter

1. hard material only