






Spring 2018

# App Features

	 <b>myControl App</b>	 <b>touchControl App</b>	 <b>easyTek App<sup>1</sup></b>
<b>Compatibility</b>	iOS/Android	iOS/Android	iOS/Android
<b>Use case</b>	Nx and Made for iPhone® hearing instruments	All primax™ Star	easyTek required
<b>Remote Control Features</b>			
Transmission signal	Bluetooth® and/or High frequency coded signal	High frequency coded signal	Bluetooth® and/or High frequency coded signal
Volume adjustment	●	●	●
Independent volume adjustment (R/L)	● <sup>1</sup>	-	-
Program change	●	●	●
Mute	●	●	●
Tinnitus signal adjustment	●	●	●
Independent tinnitus signal adjustment (R/L)	●	-	-
SoundBalance™ adjustment	●	●	●
Spatial Configurator adjustment	●	●	●
Battery status read-out	●	● <sup>2</sup>	●
Sound Exposure	●	-	-
iPhone motion sensor technology	●	-	-
Hearing Aid Airplane Mode	●	-	-
Adjustable System Sounds	●	-	-

(1) Configurable via Connexx  
 (2) Indicated by tones played in hearing instruments

## myControl™ App


The world's first use of iPhone motion sensors.

- Compatible with Pure 13 Nx, Pure 312 Nx, and Motion 13 Nx.
- Sound Exposure meter
- Stream phone calls directly into hearing instruments
- Stream music and TV directly into hearing instruments
- Adjust tinnitus therapy signal



## touchControl™ App

- Compatible with all primax devices
- Adjust Spatial Configurator
- Change programs
- Adjust volume
- Adjust SoundBalance control
- Adjust tinnitus therapy signal



## easyTek™ App

Available for Apple® iOS® and Android™ smartphones offers discreet and visual control of the hearing instruments, easyTek, and audio sources.

- Visual display: Control volume, adjust balance, select programs and audio sources, status information
- Quick access to audio sources: Stream and toggle from TV and phone



# Apps and Accessories



	myControl App	touchControl App	easyTek App	Cellion Portable Power Supply	StreamLine Mic	StreamLine TV	TV Transmitter	VoiceLink™	miniPocket™	easyPocket™	easyTek	Cellion charger	eCharger™
Pure® Charge&Go													
Pure 13 Nx													
Pure 312 Nx													
Cellion™ px	-	•	•	•		-	•	•	•	•	•	•	-
Motion® 13 Nx													
Motion 312 Nx													
Motion 13 P Nx													
Motion SX/SA px	-	•	•	-		-	•	•	•	•	•	-	• <sup>1</sup>
Motion SP px	-	•	•	-		-	•	•	•	•	•	-	-
Insio™ IIC/CIC Nx													
Insio 7px	-	•	•	-		-	•	•	•	•	•	-	-

"Made for iPhone" means that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone may affect wireless performance. The Bluetooth word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners. Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc. iPhone is a trademark of Apple Inc., registered in the U.S. and other countries. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license. Copyright © 2018 Signia GmbH. Sivantos, Inc. is a Trademark Licensee of Siemens AG. All rights reserved. 1/18 SI/18325-18