

GOVERNMENT SERVICES CUSTOM PRODUCT ORDER FORM



Life sounds brilliant.

TEL: (800) 955-4327 FAX: (732) 562-6653 EMAIL: va.customercare@sivantos.com

Account # Date:

Expedite 24-Hour Service no charge

Warranty 36-Month no charge

Bill To Address: Contact: Ship To Address:

Tel:

Email: P.O. No.

Client Information Last Name, First Name, Middle Initial

Shipping: No charge UPS 2nd Day Last 4 S.S.#

Step 1 - Side Preferences

Side Selection	Matrix Selection	
<input type="checkbox"/> Binaural	<input type="checkbox"/> Factory Select	<input type="checkbox"/> Requested Matrix
<input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> Max Matrix	L <input type="text" value="/ /"/> R <input type="text" value="/ /"/>

Hearing Thresholds

	250 Hz	500 Hz	1 kHz	2 kHz	3 kHz	4 kHz	6 kHz	8 kHz
L	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
R	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

STEP 2 - Technology and Shell Size (check selection)

	Mic	ITE (312 battery)	HS (312 battery)	ITC (312 battery)	CIC (10 battery)	IIC (10 battery)
Insio™ 7Nx (Bluetooth®)	Dir	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/>
Insio 7Nx (non-Bluetooth)	Dir	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

CROS Pure® Charge&Go Nx (not compatible with CIC or IIC)
 CROS Pure® 312 Nx (not compatible with CIC or IIC)

¹With integrated memory button in brown and mocha faceplate only

STEP 3 - Faceplate Color

	ITE	ITC/HS	CIC	IIC
Beige	<input type="checkbox"/>	<input type="checkbox"/>		
Tan	<input type="checkbox"/>	<input type="checkbox"/>		
Mocha	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Brown	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dark Brown	<input type="checkbox"/>	<input type="checkbox"/>		
Black				<input type="checkbox"/>

STEP 4 - Shell Color

	ITE	ITC/HS	CIC	IIC
Beige	<input type="checkbox"/>	<input type="checkbox"/>		
Tan	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Mocha	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Brown	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Dark Brown	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Black	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Red/Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Clear	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Deep Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

STEP 5 - Control Options

	L	R
Volume Control (N/A for CIC, IIC) ⁴	<input type="checkbox"/>	<input type="checkbox"/>
Memory Button ^{6, 7}	<input type="checkbox"/>	<input type="checkbox"/>
Telecoil ^{5, 8}	<input type="checkbox"/>	<input type="checkbox"/>
TwinPhone ⁹	<input type="checkbox"/>	<input type="checkbox"/>

⁴ only available for non-BT 7Nx
⁵ only available for Insio non-BT 7Nx
⁶ Not available for IIC
⁷ Standard for Nx BT
⁸ Not available for CIC, IIC, Insio Nx
⁹ Standard for Insio Nx

STEP 6 - Shell Options

	L	R
Retainer Loop ¹²	<input type="checkbox"/>	<input type="checkbox"/>
Removal String	<input type="checkbox"/>	<input type="checkbox"/>
Raised Volume Control Cap ¹³	<input type="checkbox"/>	<input type="checkbox"/>
Dip Canal in Soft Material ¹⁵	<input type="checkbox"/>	<input type="checkbox"/>
Canal Grip™ ¹⁶	<input type="checkbox"/>	<input type="checkbox"/>
Digital Canal Lock (ITC, CIC) ¹⁷	<input type="checkbox"/>	<input type="checkbox"/>
QuickGuard¹² (default)	<input type="checkbox"/>	<input type="checkbox"/>
Extended Receiver Tubing ^{12, 18}	<input type="checkbox"/>	<input type="checkbox"/>
Bell Bore ¹²	<input type="checkbox"/>	<input type="checkbox"/>
Clear Coat	<input type="checkbox"/>	<input type="checkbox"/>
Wax spring ¹⁹	<input type="checkbox"/>	<input type="checkbox"/>

¹² Not available for IIC
¹³ Not available for CIC & IIC
¹⁵ Not available for red, blue, clear shell color
¹⁶ Not available for clear coat, dip canal in soft material
¹⁷ Only available for HS, ITC, CIC - requires impression
¹⁸ Not compatible with WX
¹⁹ Not available with extended receiver tubing

STEP 8 - Vents (see back for vent descriptions)

	L	R
Large SAV (3.0mm) ^{10, 11}	<input type="checkbox"/>	<input type="checkbox"/>
Medium SAV (2.5mm) ^{10, 11}	<input type="checkbox"/>	<input type="checkbox"/>
Small SAV (1.6mm)	<input type="checkbox"/>	<input type="checkbox"/>
Largest Vent Possible ^{10, 11}	<input type="checkbox"/>	<input type="checkbox"/>
IROS	<input type="checkbox"/>	<input type="checkbox"/>
Semi-IROS	<input type="checkbox"/>	<input type="checkbox"/>
Pressure Vent (1.3mm)	<input type="checkbox"/>	<input type="checkbox"/>
No Vent	<input type="checkbox"/>	<input type="checkbox"/>
Micro Vent (1.0mm)	<input type="checkbox"/>	<input type="checkbox"/>

¹⁰ Not available for CIC and IIC ¹¹ Not available for gain > 55 dB

STEP 11 - CROS Pure Color Options

	CROS Pure Charge&Go		CROS Pure 312	
	L	R	L	R
Beige (1)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Granite (2)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Grey (3)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Black (5)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Silver (6)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Pearl White (7)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Golden Blonde (30)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Sandy Brown (31)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dark Champagne (37)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Dark Granite (66)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Deep Brown (72)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

STEP 9 - Canal Length


	L	R
Deep	<input type="checkbox"/>	<input type="checkbox"/>
Long	<input type="checkbox"/>	<input type="checkbox"/>
Medium	<input type="checkbox"/>	<input type="checkbox"/>
Short	<input type="checkbox"/>	<input type="checkbox"/>

(SIDE VIEW)

STEP 7 - Tragus Angle

Please determined which image best represents the wearers tragus for best microphone alignment.

Tilted tragus Vertical tragus



STEP 10 - Remote Controls

- miniPocket™
- Smart Remote

STEP 12 - Wireless Adaptors

- StreamLine Mic²¹
- StreamLine TV²¹

²¹ Nx only


If necessary, may we change the following:

Yes	No	Call
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> Faceplate (battery size)

Special instructions on back


SIGNIA CUSTOM VENT SELECTION GUIDE

IROS




- Allows low frequency sounds to enter the ear canal unaided while the high frequencies are amplified
- D-Shaped vent is stepped (or opened) at both the faceplate and canal ends
- Enables vent size to be maintained when patient has a narrow ear canal and shallow concha bowl
- No plug

Semi-IROS




- Enables vent size to be maintained when patient has a narrow ear canal
- No plug

SAV



- Diameter at the faceplate end is the same as the diameter at the canal end
 - SAV allows vent diameter to be altered using vent plugs
- Size(mm): Large-3.0, Medium-2.5, Small-1.6

No Vent



- Hearing aid does not have any venting

Special Instructions:

For internal use only

Aids: 1 2 Earmolds: 1 2 Impression: Yellow Blue Pink Peach White Purple Trans Brown Turq Green Org
Receiver: 1 2 Impressions: 1 2
Charger: 1 2 Remotes: 1 2 RECEIVED = NJ Bag # _____ ID # _____